Manual de uso / User Manual

Powerlift 135 mini-150-175

tecnimoem
care

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1.- BENEFITS, WARRANTY, TECHNICAL SERVICE

The product is intended by the manufacturer to be used in humans for the treatment or alleviation of a disease, as well as the treatment, alleviation or compensation of an injury or deficiency, according to current European regulations on **medical devices** .

The product can only be used for lifting, lowering and transferring patients.

Our products are guaranteed for a period of 3 years against any manufacturing defect (12 months in the case of batteries). In case of interventions by unauthorized personnel or failure to correctly follow the assembly, use and maintenance instructions described in this manual, all warranty and liability rights will expire. For reasons of hygiene, returns of previously used harnesses will not be accepted.

Tecmoem CARE, SL has a network of distributors for the marketing of all its products. In the event of any problem, the end user must contact the DISTRIBUTOR where the product was purchased. The distributor will have sufficient training and spare parts provided by Tecmoem to carry out the maintenance of the products.

TECNIMOEM CARE, SL defends a permanent policy of improving its products. For this reason, the characteristics of the devices can be modified without prior notice.

2.- SAFETY WARNINGS

Any person using/assembling these items must first read and understand the instructions for assembly, use and maintenance. Failure to follow these instructions may cause property damage or injury.

2.1.- Recommendations

THE LIFT SHOULD NEVER BE OPERATED BY THE PATIENT HIMSELF, but rather by a trained CAREGIVER. Familiarize yourself with the controls and safety devices before using it on patients.

Tecmoem recommends the use of the harness supplied with the crane, which has been previously tested. Otherwise, you must check that it bears the CE marking and is of the appropriate size and characteristics for the patient.

Do not use a worn harness and follow its instructions for use.

Check that the weight of the patient together with the weight of the harness does not exceed the maximum load marked on the lift (see technical specifications on page 7).

Use the loaded crane on smooth, horizontal and obstacle-free terrain.

If a patient is transferred on an inclined surface, it must not exceed a 5° slope and the help of another person is recommended.

Never leave the lift with a patient on an inclined surface.

Transfer the patient with the legs of the lift in parallel and at a speed of less than 3 km/h. (0.8m/s).

Keep the equipment away from water, humidity, splashes or corrosive environments (indoor swimming pools, the sea...). Do not use the electric lift in the shower or recharge it in the bathroom.

Do not force controls or safety devices.

Push the lift holding the front handle, never by the mast, arm, motor or the patient.

Do not raise, transfer or lower the patient with the rear wheel brakes activated.

Check that if the crane has been exposed to the sun for a long time, there may be areas that present a risk of burns.

BEFORE TRANSFERRING THE PATIENT, assess individual risks including:

- Patient mobility and holding capacity (assess type of harness for amputees, head restraint...).
- Entrapment hazard warning.
- Possible falls.
- Disturbed patients or with learning difficulties.
- Young children or pregnant women.
- People without the mental capacity to recognize unsafe actions.
- Unauthorized persons.

All persons authorized to use the lift must be able to do so in a safe and controlled manner.

The user and/or caregiver must ensure that there are no body parts in areas where there may be a risk of entrapment.

2.2.- Electrical safety precautions

This product complies with EN 60601-1 and EN 60601-1-2 standards, related to electrical devices and electromagnetic interference of medical equipment, so it does not interfere or is not interfered with when combined with other medical devices that also comply with electromagnetic regulations.

Some devices, especially older ones, that do not comply with the electromagnetic compatibility regulations may nevertheless suffer interference or cause it themselves when used with this lift. The user of such devices must ensure that possible malfunctions do not harm the patient or any other person.

Before moving the crane, make sure that the power cable is unplugged from the mains and the wheels are not stopped. Only qualified and authorized personnel can manipulate the electrical parts.

3.- ASSEMBLY

Although the use of tools is not necessary for assembly, some parts are of considerable weight. It is recommended the help of another person to handle them. Take special care to avoid possible entrapment or impact.

3.1.- Components (see page 3-5)

BASIC SET

- 1.- Frame
- 2.- Pedals
- 3.- Rear wheels (with brake)
4 legs
- 5.- Front wheels (without brake)

MAST-ARM ASSEMBLY

- 6.- Mast
- 7.- Arm

- 8.- Care handle
- 9.- Hanger
- 10.- Hooks for 4-point harness
- 11.- Electric motor
- 12.- Battery charging base (control box)
- 13.- Removable battery
- 14.- Manual control knob
- 15.- Handle (initially placed on the frame)
- 16.- Mast flange (lower engine coupling)
- 17.- Arm flange (upper coupling of the engine)

3.2.- Assembly of sets (see images on page 6)

Carefully remove the packaging.

BE CAREFUL NOT TO REMOVE THE GRAY PLASTIC SUPPLEMENTS (A)

Place the base assembly on the ground with the rear wheels stopped. Remove the handle (B) and insert the mast-arm assembly into the central frame. Finally, adjust the assembly using the handle (B), making sure that it is firm and tight. Check that the legs open and close normally and that the wheels turn correctly. Check that the hanger rotates and swings freely and is firmly attached to the arm.

The motor is supplied assembled from the factory. In its lower part it is adjusted with a screw and nut and in its upper part with a pin (C) and a safety ring (D). Make sure motor pin (J) and hand control pin (K) are FITTED FIRMLY into charging base.

4.- HANDLING INSTRUCTIONS (see images on page 6)

4.1.- Handling of the hoist

To open the legs you have to step on the left pedal. To close the legs you have to press the right pedal. To transfer patients and/or pass between doors or through narrow corridors, it is recommended that the legs are closed.

Use the control to raise or lower the crane arm. You can also use the buttons (G) on the charging base (control box). During the lowering or lifting of the patient, it is recommended that the wheels are not braked. In this way, the patient-crane unit will seek the center of gravity, avoiding discomfort to the patient.

When the arm reaches its lower or upper limit it stops automatically.

It is advisable to avoid repeatedly pressing the buttons (G), since the successive starts of the motor consume more battery, reducing its useful life.

The equipment enters standby mode if no buttons are pressed. To reactivate it, press any button.

The electrical equipment has a **red emergency stop button (F)**. With the button released the battery is ready to work. With the button pressed the battery does not act. **To release the button it must be turned clockwise and released.**

In case of battery failure, the motor has an **emergency manual lowering system**. To do this, the red wheel (E) of the motor must be pressed down, while turning clockwise to slowly lower the patient.

4.2.- Battery recharge

During lifting or lowering, the charging status lights (H) on both the charging base and the control light up green. Two green pilots indicate capacity between 66-100%, 1 green pilot indicates capacity between 33-66% and 1 orange pilot indicates capacity less than 33%.

The system emits an acoustic warning when the battery is about to run out, allowing several complete cycles to be carried out.

If you only have a removable battery and a charging base, you must bring the lift closer to a power outlet and connect the cable (L) to the charging base and then connect to the mains (110-220V).

While the battery is being recharged, the crane motor does not work, the emergency stop must not be pressed and the green light comes on. Battery stripes flash. When it is fully charged the 3 stripes will remain fixed. The green light will remain on until the network cable is disconnected.

Do not disconnect the plug from the charger by pulling on the cable as it could be damaged.

If you use the lift daily, it is recommended to recharge the battery at night. The charger does not allow the battery to overcharge. If you are not going to use the lift for a long time, it is advisable to recharge the battery at least every three months so that it is in its optimal state of use.

If you have an **additional charging base** (on request) you can recharge the battery by removing it from the charging base installed on the lift and placing it in an additional charging base (M) installed next to a power outlet.

4.3.- Battery change

To **remove the battery**, hold it by the upper handle (I) and release the metal guide from the control box support. To place it back on the charging base, rest the battery on the charging base and insert the guide again until you hear a "click".

4.4.- Overload

The control box stops working if the actuator is overloaded (due to excess load or excessive continuous operating time). It will work again after a few minutes when the overload has been removed. Approx. duty cycle: 2 minutes in continuous use and 18 minutes in rest.

4.5.- Improper uses

Improper use of this lift may result in injury to users and/or damage to equipment, as well as void product warranty. Examples:

- Use by a person who has not read this manual or who has not been trained by an authorized team member.
- Use of the crane, of the functions, of the accessories or of the movement of the crane by people not trained to use it safely.
- Use of electrical functions by several people at the same time.
- Use of the crane with a load greater than that identified on the crane.
- Irregular connection to the electrical network when charging the battery.
- Connection of devices to the crane not authorized by the manufacturer.
- Displacement of the crane pulling the power cable.

- Washing with too much water, with jet pressure or in a washing tunnel.
- Outdoor use or to transport a patient in a vehicle.
- Displacement on soft unconditioned ground.
- Movement with a patient on terrain with an inclination greater than 5°.
- Extreme use of the motor without respecting the maximum load indicated on its label.
- Use of accessories or equipment other than those recommended by the manufacturer.
- Any other use of the crane not in accordance with its intended purpose.
- Failure to respect the work cycle of the equipment: 2 minutes in continuous use and 18 minutes at rest.

5.- SLING (see images on page 6 and 7)

Tecnimoem recommends the harness supplied with the lift as it can be used by physically dependent patients, although they must have a good degree of head control. It allows easy access to technical and auxiliary personnel, and can be used both sitting and lying down.

100% Polyester Composition. The materials used to manufacture the harness are hypoallergenic . No skin reactions have been described, however, if any, it is recommended to consult the prescriber. Adjustment in three heights with safety anchor (N).

5.1.- Sling safety and care

The choice of the appropriate size as well as the correct placement of the product in the area to be treated are essential to effectively fulfill its function, for its useful life and for its safe use.

SECURITY. Never use with loads higher than the maximum load of the crane (marked on it).

Check that the sling is in perfect condition before using it on patients. It is recommended to use it previously in healthy people to familiarize yourself with it and check that it adapts to the size and characteristics of the patient.

Do not use the harness if it shows signs of tearing, wear or deformation.


The suspension straps at the ends of the harness straps (blue color) must not be used. They are safety ties in case of breakage of any of the other suspension straps (N).

CARE. Recommended washing at 30°C with neutral soap. Do not bleach or dry clean. Does not support ironing. Dry on flat surface. Do not expose to the sun or heat source. Does not support dryer.

5.2.- Placing the sling on a seated patient

Lean the patient forward, ensuring that they are well supported to reduce the risk of falling and to convey safety.

Slide the sling under the patient's back, in the center, with the opening symmetrical to the axis of the spine.

Rest the patient's back against the sling (graphic 1). Pass the leg support to the front of the chair on each side of the thigh and then pass each leg support under each thigh (chart 2)

5.3.- Placing the sling lying on a bed or on the floor


If the patient can get up from the bed or the floor, this will facilitate use. Lean the patient forward making sure they are well supported to reduce the risk of falling and to promote a sense of security.

Slide the sling under the patient's back, in the center, with the opening symmetrical to the axis of the spine.

Rest the patient's back against the sling (graphic 3).

If the patient does not have sufficient trunk control, roll up the sling halfway, turn the patient on their side and slide the sling off the raised side. Make sure the center of the harness is symmetrically located along the spine. Turn the patient the other way, unroll the sling, and return the patient to your position on it (graphic 4).


Bend your knees, one at a time, to place the harness under your legs.

5.4.- Using the leg straps

Option 1.- Leg straps crossing: For maximum security, cross the leg straps. This will reduce the risk of pitching forward and will also maintain better alignment of the hips and lower limbs (graph 5).

Option 2.- Place the leg straps together: If the patient is an amputee or feels more comfortable, pass each strap individually under both legs. Make sure the patient is not prone to spasm or there is a risk of pitching forward (Chart 6).

Option 3.- Separate leg straps: Apply each support under each leg, without crossing straps. This may be recommended when the patient has abdominal tenderness or requires intimate hygiene (graph 7).

5.5.- Fastening the suspension strips to the lateral hooks of the hanger

Fasten the shoulder and head straps before the leg straps. The harness can be used on hangers with 2 or 4 attachment points. For safety and comfort reasons the width of the hanger should correspond to the size of the harness.

The position of the patient in the harness can be adjusted by using different straps. Short shoulder straps and long leg straps produce a more upright lift, which will help with transfer to a chair. By lowering the shoulder straps and shortening the leg straps, you can achieve a reclined position, more useful for transferring from chair to bed.

6.- CLEANING AND MAINTENANCE OF THE CRANE

6.1.- Cleaning and disinfection

- Make sure that the crane is immobilized and disconnected from the electrical network and with the emergency stop pressed.
- Never wash the crane with abundant water, or with a high pressure jet, or in a washing tunnel.
- Do not use water at a temperature higher than 60° C.
- Avoid any excess water on the connection pins.
- Carefully dry the lift before using it again.
- Use clean single-use or recyclable rags.
- Use a diluted detergent or disinfectant solution or a disinfectant spray.

6.2.- Safety standards for maintenance

Maintenance operations may only be carried out by qualified and authorized technical personnel. Contact your local Tecmoem Care, SL distributor for any specific maintenance and/or technical service problem.

Before carrying out any repair or maintenance intervention:

- Make sure that the crane is immobilized (if no movement is planned) and disconnected from the electrical network.
- Block the electrical functions by pressing the red stop button.
- Never open or drill an electric motor.

6.3.- Preventive maintenance

In addition to the recommendations described in chapter 2, it is advisable to carry out an annual review in order to keep the crane and its accessories in good condition. Especially:

- Mechanisms and drive controls.
- The joints of the different movements of the crane and its accessories.
- The condition of electrical cables and the tightness of electrical appliances.
- The state of the harness (possible wear, breakage, deformation...).

It is recommended to periodically grease the joints of the mobile parts of the crane.

The periodicity and frequency of the technical revisions must be adapted to the state of the crane and to the particular conditions of use.

6.4.- List of spare parts

On request, the electrical equipment (battery, motor, controller, charging cable or control box), harness, wheels and other components of the crane's metal structure can be replaced.

6.5.- Transport and storage

During the transfer of the crane, it must be in the low position, with functions deactivated and wheels blocked and secured with straps, protected against the ingress of water or humidity (< 75%) and at a temperature between -20° C and +50° C.

7.- PROBLEM SOLVING

ISSUE	CAUSE	SOLUTION
the engine does not work	Red emergency stop button pressed	Unlock button by turning clockwise
	The battery is charging to the mains	Unplug black battery recharge cable
	Low battery	Recharge the battery
	Battery not making contact	Remove the battery and put it back
	Poorly connected control or motor cables	Connect wires correctly all the way
	Controller or motor not working	Replacing the drive or motor
Batteries not charging	Power failure	Check mains power
	Malfunction of the battery or its charging base	Replacing the battery or charging base
Engine stop	Overweight on the crane	Reduce load. Wait a few minutes
	Low battery	Recharge the battery
noises in the joints	Need for lubrication	Joint lubrication
The mast has a lot of clearance with respect to the frame	Gray plastic parts have been removed or are broken	Fit/replace the gray plastic parts on the frame
The wheels do not turn or brake properly	Dirt	wheel cleaning

tecnimoem care

SNR: ES-MF-000019199
(EU) 2017/745 MDR

TECNIMOEM CARE, S.L.
Ctra. Logroño - Mendavia, Km. 5
31230 Viana - Navarra – España (Spain)
T.: (+34) 948 646 213
www.tecnimoem.com